




OPTIMA® GC columns



tested columns for a reliable analysis

High performance OPTIMA® GC-columns, now with integrated pre-columns

For a noticeable longer lifetime, even for contaminated or matrix containing samples

- ✓ longer lifetime of the main column
- ✓ separation efficiency is not decreased by cutting or shortening the capillary
- > **saves money**
- ✓ samples with critical matrix do not endanger the main column
- > **increases flexibility**
- ✓ tight and inert connection is guaranteed, even at transport, storage, installation, dismounting, etc.
- > **increases safety of analysis**
- ✓ installation of guard columns must no longer be done personally
= "ready to use" product
- > **saves time**

Program: all main columns with 10 m guard column (with fitting deactivation)*

Cat.-No.	description	
725820.30V1	OPTIMA-5 MS-Accent	0,25 mm ID, 0,25 µm film, 30 m
726220.30V1	OPTIMA-5 MS	0,25 mm ID, 0,25 µm film, 30 m
726056.30V1	OPTIMA-5	0,25 mm ID, 0,25 µm film, 30 m
726600.30V1	OPTIMA WAX	0,25 mm ID, 0,25 µm film, 30 m
725820.60V1	OPTIMA-5 MS-Accent	0,25 mm ID, 0,25 µm film, 60 m
726220.60V1	OPTIMA-5 MS	0,25 mm ID, 0,25 µm film, 60 m
726056.60V1	OPTIMA-5	0,25 mm ID, 0,25 µm film, 60 m
726600.60V1	OPTIMA WAX	0,25 mm ID, 0,25 µm film, 60 m

* further combinations on request

www.mn-net.com

MACHERY-NAGEL

MACHERY-NAGEL GMBH & CO. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany

Germany

and international:

Tel.: +49 (0) 24 21 96 90

Fax: +49 (0) 24 21 96 91 99

e-mail: sales-de@mn-net.com

Switzerland:

MACHERY-NAGEL AG

Tel.: +41 (0) 62 388 55 00

Fax: +41 (0) 62 388 55 05

e-mail: sales-ch@mn-net.com

France:

MACHERY-NAGEL EURL

Tel.: +33 (0) 3 88 68 22 68

Fax: +33 (0) 3 88 51 76 88

e-mail: sales-fr@mn-net.com

USA:

MACHERY-NAGEL Inc.

Tel.: +1 484 821 09 84

Fax: +1 484 821 12 72

e-mail: sales-us@mn-net.com

